**heh**<sup>®</sup>  
instruments

**TechPatient**<sup>®</sup>  
**CARDIO**

**VERSION 3.2**

## **Operator's Manual**

*The described device is to be used by technically skilled personnel.*

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## FCC Exemption / Exención FCC

The device is exempt from FCC regulations under Part 15, Section 103 c.

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## Electromagnetic compliance information Información de compatibilidad EM

While this simulator is compliant with the applicable regulatory EMC requirements, there is no guarantee that interference will not occur in a particular installation. To minimize the potential for the hardware to cause interference to radio and television reception or to experience unacceptable performance degradation, install and use this hardware in strict accordance with the instructions in this operator's manual.

If this hardware does cause interference with licensed radio communications services or other nearby electronics, which can be determined by turning the hardware off and on, you are encouraged to try to correct the interference by one or more of the following measures:

- Reorient the antenna of the receiver (the device suffering interference).
- Relocate the transmitter (the device generating interference) with respect to the receiver.
- Plug the transmitter into a different outlet so that the transmitter and the receiver are on different branch circuits.

Some hardware may require the use of a metal, shielded enclosure (windowless version) to meet the EMC requirements for special EMC environments such as, for marine use or in heavy industrial areas. Refer to the hardware's user documentation and the DoC for product installation requirements.

When the hardware is connected to a test object or to test leads, the system may become more sensitive to disturbances or may cause interference in the local electromagnetic environment.

Operation of this hardware in a residential area is likely to cause harmful interference. Users are required to correct the interference at their own expense or cease operation of the hardware.

Changes or modifications not expressly approved by HE Instruments could void the operator's right to operate the simulator under the local regulatory rules.

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A patient simulator is an electronic device that through the simulation of specific parameters in preset conditions, is used for testing, verification and calibration of diagnostic equipment.

The TechPatient CARDIO simulator is specifically designed to test 12-leads ECG diagnostic equipment, such as ECG machines, patient monitors and Holter recorders among others.

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## The 12-Lead ECG El ECG de 12 Derivaciones

The standard 12-lead electrocardiogram is a representation of the heart's electrical activity recorded from ten electrodes on the body surface. It provides spatial information about the heart's internal activity, each of the 12 leads representing a particular orientation in space.

It's an improvement from the original 3 and 5-leads ECG, therefore, 12-leads simulation is compatible with 3 and 5-leads diagnostic equipment.

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The TechPatient CARDIO simulator connects to diagnostic equipment through its 10 biopotential multi-connectors. Each connector is compatible with 4mm snap, 4mm banana and 2mm pin patient cables. To connect banana patient cables, unscrew the metallic cap to expose the 4 mm jack.

A hi-signal alligator cable is also included. It delivers a 200x preamplified signal for premature circuit testing during electronic design. This output is DC coupled to ease connection to analog-digital converters.

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
The TechPatient CARDIO simulator is a low voltage battery powered device that presents no electrical safety risk to its operator. However, when connected to diagnostic equipment in fault condition, care should be taken to avoid electrical shocks conducted through the patient cable.

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





Pulse waveform at 60 PPM  
Onda de pulso a 60 BPM


Sinusoidal waveform at  
Onda senoidal


Triangular waveform at 2 Hz  
Onda triangular a 2 Hz

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

<sup>1</sup> Calibrated accuracies only apply to calibrated devices for the 12 months following the purchase of a calibration certificate.

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Hold for more than 2 seconds to turn the device on and off.

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When browsing parameters, it enters the current to change its value.

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The up and down keys navigate the list of parameters while browsing or edit the current value. The pointer indicates the current position.

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## Menu options / Opciones del menú

**Mode:** **ECG:** Generates 12-Lead realistic cardiac waveforms.

**Sin:** Generates three-phase sinusoidal waveforms on 3 leads and single-phase sinusoidal waveforms on the precordial leads.

**Tri / Quad:** Generates two-phase triangular / square waveforms referenced to a common lead.

**Pulses:** Generates a two-phase pulsed waveform with the same features as the cardiac waveform. Waveforms are referenced to a common lead.

**Base:** Generates a baseline of a predefined mean value.

**Config:** Enables the selection of interface language and line frequency.

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: Indicates the occurrence of a heartbeat. It flashes a ! symbol when heartbeats are chaotic.

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**Frequency:** Sets the waveform Frequency in 5 Hz increments.

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**Amplitude:** Sets the waveform mean Amplitude in mV (referenced to Lead II).

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**% Base:** Sets the percentage of baseline.

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**Noise:** **No:** Leaves the generated waveform unchanged.

**Line:** Adds line noise to the waveform (50 Hz / 60 Hz).

**Ripple:** Adds ripple noise to the waveform (100 Hz / 120 Hz and harmonics).

**Random:** Adds random noise with a near gaussian distribution.

**BLWand:** Adds baseline wander, usually associated with patient's motion or ventilation.

**Tremor:** Adds noise associated with muscle tremor.

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**Temblo:** Agrega ruido asociado al temblor muscular.

**Common:** Selects the reference common lead when generating two-phased waveforms.

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**Light:** Select whether the lights must remain on, off or start only when touching the keypad.

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**Sleep:** Selects whether the device must switch off after 10 minutes without user intervention.

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**Battery:** Indicates the remaining percentage of battery charge.

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The following explains a generic verification procedure for ECG based diagnostic equipment. Please refer to your equipment's service manual for a detailed and specific procedure.

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HE Instruments warrants the TechPatient CARDIO simulator from manufacturing defects in materials or workmanship, when properly maintained, and under normal use and service, which means the measurement of any of its outputs without applying external signals to them. The application of external signals such as a defibrillator shockwave is explicitly prohibited and will expire this warranty.

HE Instruments' sole obligation under this warranty shall be limited to the repair or replacement, at its option, of any defective part of said Patient Simulator which is the result of defective material or defective workmanship of parts installed by HE Instruments.

This warranty will expire 6 months from the date of delivery, and repairs under this warranty shall be at repair facilities designated by HE Instruments. Transportation expenses to and from the repair facility are to be born by the purchaser.

The warranty herein is expressly in lieu of any and all other warranties expresses or implied. No implied warranty of merchantability is made and there are no warranties which extend the use of the Patient Simulator for purposes beyond the description hereof.

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To replace the battery, open its rear door, remove the battery from its holder, and replace it with a new one. The device should never be used while the battery door remains open. Only replace it with a PP3 compatible 9V battery. Dispose discharged batteries in the appropriate environmentally clean way, according to your government regulations. Moisture should only be removed with a piece of dry cloth. Avoid chemicals and other wet cleaning products.

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---

---

For any further information regarding this or any other HE Instruments product contact us at:

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*All rights reserved © 2011, HE Instruments*  
*Todos los derechos reservados © 2011, HE Instruments*